For in vitro diagnostic use only. For professional use only.

Cat No: BS-STI-T-25 / BS-STI-T-100

Sexually Transmitted Infections RT-qPCR Panel



Package Insert

Component	Int	25 Reactions	100 Reactions	
2X Prime Script Mix	Optimized ready-to-use mix for RT-qPCR assay		1 x 1000 μL	4 x 1000 μL
STI Oligo Mix 1	Specific nucleic acid amplification and detection: FAM: Herpes Simplex Virus-1 HEX: Mycoplasma genitalium CY5: Herpes Simplex Virus-2		1 x 125 μL	1 x 500 μL
STI Oligo Mix 2	HEX: Human (Internal Control-IC) ROX: Neisseria gonorrhoeae		1 x 125 μL	1 x 500 μL
STI Oligo Mix 3	FAM: Treponema pallidum HEX: Trichomonas vaginalis ROX: Ureaplasma parvum/ urealyticum		1 x 125 μL	1 x 500 μL
STI Oligo Mix 4	FAM: Mycoplasma hominis HEX: Haemophilus ducreyi ROX: Streptococcus agalactiae CY5: Gardnerella vaginalis		1 x 125 μL	1 x 500 μL
5TI 1 / PC-STI 2 / PC-STI 3 / PC-STI 4	Positive Control (PC)		1 x 100 μL	1 x 100 μL
NTC	Negative Control (NTC)		1 x 1000 μL	1 x 1000 μL
e 2. Transport Condition, Storage Condi	tion, and Shelf Life of the Components			
Component	Transport Condition	Storage Condition*		Shelf Life
2X Prime Script Mix		(-22) °C – (-18) °C	(-22) °C – (-18) °C	
Oligo Mix		(-22) °C – (-18) °C	(-22) °C – (-18) °C	
PC	(-22) C - (+8) C	(-22) °C – (-18) °C before opening, (+2) °C – (C before opening, (+2) °C – (+8) °C after first thaw	
NTC	(-22) °C – (-18) °C before opening, (+2) °C – (-		+8) °C after first thaw	

Table 3. Components Required but Not Included with The Test

Components Required but Not Included with The Test

1. Magnetic Induction Cycler (Mic) (Bio Molecular System - BMS) or/and CFX96 Touch™/CFX96™ Dx/CFX Opus 96™/CFX Opus 96™ Dx (Bio-Rad) Real-Time PCR systems

 $\textbf{2.} \qquad \text{Micropipettes and compatible filtered pipette tips (nuclease-free) suitable for transferring 1-10, 10-100, and 100-1000 \, \mu L of liquid liquid$

3. A centrifuge or Mini-spin

4. Vortex

5. Reaction tubes and caps/films specific to RT-qPCR instruments and compatible with the reaction volume

Table 4. Intended Use, Test Principle, and Analytical Specifications				
Function	Aid to diagnosis	Sample Type(s)	Table 5	
Analyte	Table 1	Nucleic Acid Preparation Method(s)	Table 5	
Qualitative/Quantitative	Qualitative	Validated PCR Instruments	Table 3	
Test Principle	Reverse Transcription and Real-Time PCR (RT-qPCR)	Inclusivity and Exclusivity	Validated on the reference strains and the field isolates	
Automated/Manual	Manual			
Intended Users	Professional use	Limit of Detection (LoD)	Table 5	
Target Population	Individuals with the suspected infection	Sensitivity and Specificity	98.76% and % 99.28%	
Table 5. Collection, Storage, and Transfer of Clinical Specimens / Nucleic Acid Preparation Methods and the Respected LoD Values				

Sample Type**	Sample Transfer	Sample Storage	Nucleic Acid Preparation Method	LoD (cp/mL)
Throat swab, vaginal swab, urethral swab, urogenital swab, endocervical swab samples	vNAT® Transfer Tube (Cat. No: BS-NA-513m)	3 months at (+2) °C – (+8) °C 1 year at -20 °C	Nucleic acid preparation is not needed. The samples can be used directly in RT-qPCR.	400-800
Urine samples	Viral Transport Medium (VTM) (CDC SOP#: DSR-052-05, without antibiotics)	RINA™ M14 Nucleic Acid Extraction Device 3 days at (+2) °C - (+8) °C 1 uppert 70 °C 1 uppert 70 °C		100-200
	Preservative-free sterile tubes/containers	i year at -70°C	(Robot Model No: EXM3000, Kit Cat. No: ZFNAE01)	100-200

** Clinical specimens should be collected by a healthcare provider in accordance with national/international clinical specimen collection regulations.

For in vitro diagnostic use only. For professional use only.



1. RT-qPCR Application Protocol

Before starting the assay, please consider the following:

- 1. The kit was validated only for the template nucleic acid volume which is 25% of the total RT-qPCR volume.
- 2. The kit cannot be used with real-time PCR instruments without periodic maintenance records.
- 3. The kit for Bio-Rad Real-Time PCR systems has been validated with white reaction tubes specific to these systems. Clear reaction tubes result in 5-10 times lower fluorescence signal in Bio-Rad instruments compared to white reaction tubes. In addition, device-specific reaction tubes should be used in the BMS device. The kit's stated analytical performance can only be achieved using validated tubes.
- 4. To test for contamination, a negative control reaction containing NTC (Nuclease-free Water) must be set up in each run.

Program the RT-qPCR device as follows and add the reagents into the RT-qPCR tubes, close the tubes, place them into the RT-qPCR instrument and start the run (Table 6)

Table 6. Reaction Setup and Real-Time PCR Program RT-qPCR Program QR Code for Thermal Protocol and Plate Setup CFX96 Touch™/CFX96™ Dx/CFX Opus 96™/CFX Opus 96™ Dx (Bio-Rad) and Reaction Setup Magnetic Induction Cycler (Mic) (Bio Molecular System - BMS) Reagent Volume/Rxn Step Cycle No. Temperature Duration **Reverse Transcription** 1 Cycle 52 °C 3 min 2X Prime Script Mix 10 µL Pre-Incubation 1 Cvcle 95 °C 10 sec Denaturation 12 Touchdown 95 °C 1 sec Cvcles: 1 °C decrement in Oligo Mix Annealing and 5 uL 67 °C to 56 °C annealing 15 sec Extension temperature per cvcle Denaturation 95 °C 1 sec Template Nucleic 5 μL Annealing and Acid/NTC/PC 55 °C 15 sec Extension 30 Cycles **Total Reaction** (FAM-Green)/(HEX-Yellow) https://www.bioeksen.com.tr/files/L TD 43P/ 20 uL **Detection (Reading)** Volume (ROX-Orange)/(CY5-Red)

WARNING: The RT-qPCR thermal programs (Bio-Rad and BMS) and the plate setup (Bio-Rad) file should be downloaded from the QR code or link above.

2. Interpretation of the Assay Results

All default analysis options (e.g., auto-calculated threshold) in the PCR instrument software should not be changed to calculate Cq values.

The shape of the amplification curves should be examined for all reaction wells returning with Cq values. All the sigmoidal curves above the threshold should be recorded as "positive" and their Cq values should be recorded. Non-sigmoidal curves should be recorded as "negative".

Table 7. Expected Performance of Kit Controls

Control Type	Dumpere	Expected Results and Cq Values		
	Purpose	IC (HEX)	Target	
Negative Control	Contamination control during RT-qPCR	Not detected (No Cq)	Not detected (No Cq)	
Positive Control	Reagent integrity	Detected (Cq≤30)	Detected (Cq≤30)	
Internal/Extraction Control	To monitor the integrity of nucleic acid extraction and RT-qPCR from each sample	Detected (Cq≤30) If the IC is "Not detected", check the target Cq.	If the target is " Detected " according to the result interpretation criteria, IC is valid.	

If any control does not work as described above, the run is reported as follows:

- 1. Contamination: If Cq≤30 in any NTC test channel.
- Recommended action: Repeat the analysis paying attention to the "Warnings and Limitations" section.

Reagent Problem: In case a sigmoidal curve with a Cq≤30 cannot be obtained for any of all the samples tested in the run, including the controls.
Recommended action: If all the tested samples show negative results for the target pathogens and controls, the run is considered invalid. In this scenario, it is essential to conduct testing on the "Positive Control(s)" included in the kit. A negative Positive Control test result suggests a potential "Reagent Problem." If encountered, please reach out to the manufacturer for further assistance.

Invalid: If the IC (Internal Control) and all other targets are "Not detected".
Recommended action: Sampling isn't successfully done, or there is a problem during the sample transportation. A new sample from the same patient should be collected and tested again.

If all the controls are valid, the results are interpreted as follows:

Table 8. Interpretation of Patient Results

Target	Internal Control (IC)	Result Interpretation		
Positive (+)	Positive (+) or Negative (-)	Results are valid Target is detected	If 26 <cq "low="" positive"<br="" ≤30="">If 16<cq≤26 "positive"<br="">If Cq≤16 "High Positive"</cq≤26></cq>	
Negative (-)	Positive (+)	Results are valid Target is not detected		

The results generated by the qPCR instruments can be reported manually, as explained earlier, or automatically using the "Sigmoida" software. To obtain the "Sigmoida" software installer, please send an email to support@bioeksen.com.tr.



For in vitro diagnostic use only. For professional use only.

1. Warnings and Limitations

- 1. False-negative results may occur if inadequate numbers (lower than the LoD) of organisms are present in the specimen.
- 2. Mutations within the target regions could affect primer and/or probe binding, resulting in failure to detect the presence of agents.
- 3. The use of cotton or calcium alginate swabs or swabs with wooden sticks can lead to false negative results since they may contain substances that inactivate some viruses and inhibit PCR.
- 4. A false-negative result may occur if a specimen is improperly collected, transported, or handled.
- 5. The clinical specimens shall be collected by a healthcare provider in accordance with the specimen collection guidelines.
- 6. Test procedures should be performed by personnel trained in the use of the kit.
- 7. Except for liquid transfers, sample tubes should always be kept closed.
- 8. Filtered and nuclease-free pipette tips should be used for sample transfer.
- 9. The components in the kit should not be used together with different lot numbers or chemicals of the same name but from different manufacturers.
- 10. The caps of the reaction tubes must not be opened after the PCR run. The PCR tubes should be placed in a bag and thrown away after the bag is tightly closed.
- 11. The surfaces of the workbenches should be wiped with freshly diluted 10% bleach (0.5% NaClO) at the beginning and end of each day.
- 12. Disposal of waste must be carried out in accordance with local, state, and federal regulations.

2. Explanation of Symbol

Symbol	Title of Symbol	Symbol	Title of Symbol	Symbol	Title of Symbol
CE	European Conformity CE Mark	LOT	Batch code	×	Keep away from sunlight
IVD	In vitro diagnostic medical device	REF	Catalogue number	***	Protect from heat and radioactive sources
	Manufacturer	NON	Non-sterile		Do not use if package is damaged and consult <i>instructions for use</i>
\square	Use-by date	- 	Consult instructions for use or consult electronic instructions for use	Ť	Keep dry
CONTROL -	Negative control	\triangle	Caution	<u><u>†</u>†</u>	Keep upright
CONTROL +	Positive control	X	Temperature limit	Σ	Contains sufficient for <n> tests</n>
CONTROL	Control				

3. Manufacturer and Technical Support



Bioeksen AR GE Teknolojileri A.Ş

Address: Huzur Mah. Metin Oktay Cad. Nurol Life Sitesi D Blok No:3/31, 34396 Sarıyer/İstanbul-TÜRKİYE

Phone: +90 (212) 285 10 17, Fax: +90 (212) 285 10 18

Web: www.bioeksen.com.tr, e-mail: info@bioeksen.com.tr,

Technical Support: support@bioeksen.com.tr

Notice to User: Please inform us about product-related incidents at "vigilance@bioeksen.com.tr" within 24 hours.

ALL RIGHTS RESERVED